



A large, stylized red shape in the top left corner, resembling a partial circle or a speech bubble.

# Safety Blood Collection

Welcome to Caresfield Safety Blood Collection – where precision meets safety, and patient care is paramount. From straight and butterfly needles to safety holders, our innovative solutions are engineered to streamline blood collection while ensuring comfort and reliability. Explore how Caresfield is elevating blood collection standards.

A large, stylized blue shape on the right side, resembling a partial circle or a speech bubble.

## Safety Slide Lock Butterfly Needles



Safety needles with large butterfly wings that can be easily slid down after use, securely covering the needle

## Safety Push Button Butterfly Needles



Safety needles with a quick, one-handed needle retraction mechanism, designed for easy activation

## Safety Straight Needles



Safety needles with a hinged safety shield that covers the needle after use, featuring a large finger grip for control

## Needle Holders & Transfer Devices



Needle holders and transfer devices for safely transferring to syringes, empty tubes, or blood culture bottles, compatible with all Vacutainer tubes and blood collection needles

# caresfield®



*Real people. Real answers.*

**For over 30 years, Caresfield has delivered critical identification and specimen collection products to health care facilities across the U.S. Customers love us for our short lead times, same-day shipping and exceptional service. The Caresfield difference includes:**

- Same-Day Shipping – Orders placed before 2 PM CST ship the same day.
- Easy Online Ordering – Quick purchasing at your convenience.
- Human Customer Support – Real people, ready to assist when you need it.
- One-Stop Shop – A wide selection of supplies - all from a single trusted source.

## Safety Slide Lock Butterfly Needles

Our Safety Slide Lock Butterfly Needles are designed with phlebotomists in mind. The innovative sliding safety shield provides protection against needle injuries, while the large, ergonomic butterfly wings ensure a stable and comfortable grip. The translucent shield enhances blood flashback visibility, giving you clear confirmation of successful vein access. Experience smooth, reliable, and secure collections with every use!



SKU #61636

Item #	Description	QTY/UOM	UOM
61634-CS	21G Precision Safety Slide Lock Butterfly Blood Collection Set, 0.75in needle, 12" tube, Sterile, 50/BX - 4 BX/CS	200	CS
61635-CS	23G Precision Safety Slide Lock Butterfly Blood Collection Set, 0.75in needle, 12" tube, Sterile, 50/BX - 4 BX/CS	200	CS
61636-CS	25G Precision Safety Slide Lock Butterfly Blood Collection Set, 0.75in needle, 12" tube, Sterile, 50/BX - 4 BX/CS	200	CS
61340-CS	21G Safety Slide Lock Butterfly Blood Collection Set, 0.75in needle, 12" tube, Sterile, 50/BX - 4 BX/CS	200	CS
61341-CS	23G Safety Slide Lock Butterfly Blood Collection Set, 0.75in needle, 12" tube, Sterile, 50/BX - 4 BX/CS	200	CS
61342-CS	25G Safety Slide Lock Butterfly Blood Collection Set, 0.75in needle, 12" tube, Sterile, 50/BX - 4 BX/CS	200	CS

## Safety Push Button Butterfly Needles

These needles make blood collection safer and more comfortable. With a single-handed push, the needle retracts to reduce sharp exposure. The unique bevel needle design increases blood flow for a smoother, more comfortable patient experience.



One-hand needle retraction

SKU #61343-CS

Item #	Description	QTY/UOM	UOM
61343-CS	21G Precision Safety Push Button Butterfly Blood Collection Set, 0.75in needle, 12" tube, 50/box - 4 BX/CS	200	CS
61344-CS	23G Precision Safety Push Button Butterfly Blood Collection Set, 0.75in needle, 12" tube, 50/box - 4 BX/CS	200	CS
61345-CS	25G Precision Safety Push Button Butterfly Blood Collection Set, 0.75in needle, 12" tube, 50/box - 4 BX/CS	200	CS
61350-CS	21G Safety Push Button Butterfly Blood Collection Set, 0.75in needle, 12" tube, 50/box - 4 BX/CS	200	CS
61351-CS	23G Safety Push Button Butterfly Blood Collection Set, 0.75in needle, 12" tube, 50/box - 4 BX/CS	200	CS
61352-CS	25G Safety Push Button Butterfly Blood Collection Set, 0.75in needle, 12" tube, 50/box - 4 BX/CS	200	CS

## Safety Straight Needles

Caresfield's straight needles are designed for safety and ease of use. The hinged safety shield with an extra-large finger grip helps protect against needle injuries, while built-in needle orientation ensures precise positioning. Plus, the audible click lets you know the safety shield is locked in place—giving you confidence with every draw.



SKU #61625

Item #	Description	QTY/UOM	UOM
61625-CS	21G Green Safety Blood Collection Needle w/ integrated safety shield 1.25in Sterile 48/BX - 10 BX/CS	480	CS
61626-CS	22G Black Safety Blood Collection Needle w/ integrated safety shield 1.25in Sterile 48/BX - 10 BX/CS	480	CS
61627-CS	23G Blue Safety Blood Collection Needle w/ integrated safety shield, 1.25in Sterile 48/BX - 10 BX/CS	480	CS

## Needle Holders & Transfer Devices

Caresfield's single-use holders make specimen transfers easy and hassle-free. Compatible with all Vacutainer® tubes and blood collection needles, they provide a safe, seamless way to transfer samples from a syringe to a tube or blood culture bottle—keeping your workflow smooth and efficient.



SKU #61637

SKU #61639

SKU #61641

Item #	Description	QTY/UOM	UOM
61637	Holder for Blood Tubes, Square Edge, Clear Plastic, 250/BAG	250	BG
61639	Blood Transfer Holder Device, Red Female Luer Connection, Sterile, 200/BAG	200	BG
61641	Blood Transfer Holder Device, Blue Male Luer Connection, Sterile, 200/BAG	200	BG
61637-CS	Holder for Blood Tubes, Square Edge, Clear Plastic, 250/BAG - 4 BAGS/CS	1000	CS
61639-CS	Blood Transfer Holder Device, Red Female Luer Connection, Sterile, 200/BAG - 4 BAGS/CS	800	CS
61641-CS	Blood Transfer Holder Device, Blue Male Luer Connection, Sterile, 200/BAG - 4 BAGS/CS	800	CS

Contact Us: **952.890.7100** or [sales@caresfield.com](mailto:sales@caresfield.com) to request a catalog, conversion, or more information!

**caresfield**<sup>®</sup>  
let's talk