

# Safety Lancets

Precision. Safety. Comfort.

Caresfield's disposable safety lancets are designed with safety and comfort in mind, helping to reduce the risk of needlestick injuries while providing you with a complete range of sizes and activation methods to suit the unique needs of every patient.



## Patient Safety and Comfort

- **Needlestick Prevention:** Needles and blades automatically retract before and after use, meeting OSHA requirements and preventing cross-contamination.
- **Pain Reduction:** High-speed penetration and tri-bevel needles minimize patient pain and tissue trauma.
- **Versatility:** A complete portfolio ranging from micro-preemie heel sticks to high-flow adult fingersticks.

## Optimized Workflow



- **Consolidated Supply Chain:** Source your Blood Collection Needles, Phlebotomy Tourniquets and Lancets from a single, trusted partner.
- **Intuitive Design:** Color-coded packaging and simple activation methods (Push Button or Contact Activated) reduce training time and errors.

## Heel Incision Devices (Neonatal)

Engineered with a sweeping blade motion for high blood flow with minimal depth, reducing bone injury risk in newborns.



SKU	Description	Depth	Length	Indication
64150	Incision Device (Preemie)	0.65 mm	1.4 mm	Micro-Preemie (Low birth weight)
64151	Incision Device (Newborn)	0.85 mm	1.75 mm	Newborn Metabolic Screening
64152	Incision Device (Toddler)	1.0 mm	2.5 mm	Toddler / High Volume



## Contact Activated Safety Lancets

Simply push against the skin to fire. Ideal for rapid, consistent sampling.



SKU	Gauge	Depth	Flow	Use Case
64166	30G	1.6 mm	Low	Glucose / Pediatric Fingers
64190	28G	1.6 mm	Low	Frequent Glucose Testing
64167	28G	1.8 mm	Low	Standard Glucose
64168	26G	1.8 mm	Medium	General Diagnostics
64169	23G	1.8 mm	Medium	Cholesterol / HbA1c
64162	21G	1.8 mm	High	Coagulation (PT/INR)
64170	21G	2.2 mm	High	Calloused Skin / Manual Laborers


## Push Button Safety Lancets

Button-activated for precise placement and control.



SKU	Gauge	Depth	Flow	Use Case
64174	25G	1.4 mm	Low	Pediatric / Thin Skin
64186	25G	1.6 mm	Low	Standard Glucose
64184	23G	1.8 mm	Medium	General POCT
64165	21G	1.8 mm	High	Capillary Tubes
64193	21G	2.0 mm	High	Blood Gases
64175	18G	1.8 mm	Max	High Volume Collection

## Specialty Lancets



SKU	Description	Details
64142	Adjustable Depth Lancet	28G with 3 settings (1.4 / 1.8 / 2.2 mm)
64176	Twist Cap Safety Lancet	21G, 2.0 mm

 Single-Use Only

 Not Made With Natural Rubber Latex

8095 215th Street West, Lakeville, Minnesota  
952.890.7100 | sales@caresfield.com | caresfield.com

**caresfield**<sup>®</sup>  
let's talk