

Disposable Tourniquets for Phlebotomy

Caresfield's non-latex disposable tourniquets are designed with safety and comfort in mind, helping to reduce the risk of leave-on accidents while providing you with the widest range of size and color options to suit the unique needs of every draw.



Patient Safety and Comfort



- Bright colors reduce risk of leave-on accidents
- Single-use design enhances hygiene and helps to prevent cross-contamination/infection
- Non-latex, odor free, and powder-free

Optimized Mechanical Properties



- Superior compression, excellent stretch, and tear resistant
- Specialty tourniquets are properly sized to help assure correct amount of pressure is applied while providing maximum comfort for infant, pediatric, and bariatric patients




Standard Tourniquets for Phlebotomy

80% of infections are transmitted through touch¹, so there is an added risk of cross-infection for every product that is used to treat multiple patients. Single-use tourniquets from Caresfield are the safer, cost-effective alternative to reusable solutions. Our standard, adult tourniquets measure 1" x 18" and come in more packaging, texture, and color options than any other line today!

Bright Colors to Prevent Forgotten Tourniquets


When tourniquets are left in place for too long, you risk hemolysis and damage to the blood sample, or worse, significant injury can occur to the patient in the form of circulatory, neurological, vascular, and muscular damage². The most common cause for these "forgotten tourniquet" incidences is due to environment-related factors, such as the tourniquet being found under gowns and healthcare drapes³.



So what can you do to prevent these unfortunate mistakes?

You can make sure that your tourniquets stand out and get noticed in a sea of monochrome or muted colors typically used in healthcare! By using brightly colored tourniquets, caregivers are more likely to notice a tourniquet and with five standard color options, you can easily color-code your tourniquets by department.

Satisfy Users' Grip Preference



Smooth and textured tourniquets both provide exceptional grip and excellent strength. If you want to give your patients a softer experience, try our textured tourniquets. But if you're looking for an extra dose of grip, our smooth tourniquets were designed for you!

For those who want the best of both worlds, order our Versaquet™ Multi-Texture Tourniquet (patent pending). With one smooth side and one textured side, Versaquet can satisfy multiple end users.

Packaging Choices for Every Department

Roll Pack

Notch cut every 18 inches for quick and easy dispensing



Flat Packed

Allow tourniquets to be easily dispensed and are most commonly used for primary draw stations



Dispenser Box

Provide extra protection for tourniquet rolls plus offers notch cuts for easy dispensing



Individually Rolled

Compactly rolled and held in place with a small band, these are perfect for phlebotomists on the go



1. Centers for Disease Control and Prevention (CDC). CDC Quick Reference Guide for Public Information on Infection Control.
2. Pennsylvania Patient Safety Authority. Forgotten But Not Gone: Tourniquets Left on Patients
3. Pennsylvania Patient Safety Authority. The Forgotten Tourniquet – An Update

Specialty Tourniquets

Patients come in all sizes, so why shouldn't your tourniquets? Successful venipuncture requires that a tourniquet be tied tight enough to stop arterial blood flow, causing the collection vein to swell – so using the proper tourniquet with the right mechanical properties for a patient is absolutely necessary.

Orange Bariatric Tourniquets

For use on your largest patients, bariatric smooth tourniquets measure 1.5" x 32", nearly twice the length and 1.5X the width of standard, adult tourniquets.



No more tying



THE CHALLENGE
When standard tourniquets are not long enough for larger patients, nurses sometimes tie two tourniquets together to extend the length.
THE SOLUTION
Caresfield bariatric tourniquets are nearly twice as long as a standard tourniquet so that you can comfortably apply the band to a larger patient.
THE RESULT
Bariatric tourniquets provide you with process efficiency while respectfully providing your patients with the comfort that they deserve.

Pink Pediatric Tourniquets

Designed for most pediatric patients, from toddlers to school-aged children, our pediatric smooth tourniquets measure 0.75" x 12" and are specially designed to provide just enough grip and pressure for little ones.



No more cutting



THE CHALLENGE
Standard tourniquets may be oversized for pediatric patients, so at many facilities today, nurses might cut tourniquets to a shorter length.
THE SOLUTION
Caresfield pediatric tourniquets are sized for smaller patients, allowing you to easily tie a tourniquet with the appropriate amount of pressure for venipuncture.
THE RESULT
Pediatric tourniquets eliminate the need for cutting, speeding up the venipuncture process, and reducing the risk of cross infection from a scissors.

Infant Tourniquets

Made for your most adorable patients, these infant tourniquets measure a wee 0.5" x 9" and offer a smooth, comfortable texture. Infant tourniquets come in unique "baby"-blue or purple colors so you can easily identify the right product for your tiniest patients.



No more rubber bands




THE CHALLENGE
NICU nurses, lacking an appropriately-sized solution for infant patients, might use rubber bands or cut off glove tips for tourniquets.
THE SOLUTION
Caresfield infant tourniquets are sized appropriately for infant limbs, come in accessible packaging, and have easily identifiable colors.
THE RESULT
Infant tourniquets improve process efficiencies for NICU nurses, increase patient comfort, reduce the risk of latex allergies, and reduce the risk of cross-contamination due to non-tourniquet materials being used incorrectly.

Tourniquets for Every Patient

With more options than other brands of disposable tourniquets, including three textures, five sizes, five color options, and four packaging options, you are sure to find your perfect fit! Caresfield has a tourniquet for every type of patient you might see.

Don't see what you are looking for?

Don't worry, we can make it for you. Below are our most popular items. If you would like a different color, size, texture, or packaging option, just let us know!



Ordering Information

See below for ordering information of our standard products and contact us at sales@caresfield.com or **952.890.7100** to place an order or inquire about a custom tourniquet.

Standard/Adult 1" x 18" Non-Latex Tourniquet							
Texture	Packaging	Orange	Blue	Magenta	Green	Yellow	Assorted
Smooth	Flat Packed 1000/case	48NLO-FP	48NLB-FP	48NLM-FP	48NLG-FP	48NLY-FP	48NLA-FP
	Roll Pack 25/roll, 40 rolls/case	48NLO-RL	48NLB-RL	48NLM-RL	48NLG-RL	48NLY-RL	48NLA-RL
	Dispenser Box 25/roll, 40 boxes/case	48NLO-DBX	48NLB-DBX	48NLM-DBX	48NLG-DBX	48NLY-DBX	48NLA-DBX
	Individually Rolled 1000 rolls/case	48NLO-IR	48NLB-IR	48NLM-IR	48NLG-IR	48NLY-IR	48NLA-IR
Textured	Flat Packed 1000/case	48TXO-FP	48TXB-FP	48TXM-FP	48TXG-FP	48TXY-FP	48TXA-FP
	Roll Pack 25/roll, 40 rolls/case	48TXO-RL	48TXB-RL	48TXM-RL	48TXG-RL	48TXY-RL	48TXA-RL
	Dispenser Box 25/roll, 40 boxes/case	48TXO-DBX	48TXB-DBX	48TXM-DBX	48TXG-DBX	48TXY-DBX	48TXA-DBX
	Individually Rolled 1000/case	48TXO-IR	48TXB-IR	48TXM-IR	48TXG-IR	48TXY-IR	48TXA-IR
Versaquet (Multi-Texture)	Flat Packed 1000/case	48DSO-FP	48DSB-FP	48DSM-FP	48DSG-FP	*	*

Special Sizes							
Texture	Packaging	Infant (Blue) 0.5" x 9"	Infant (Purple) 0.5" x 9"	Infant (Assorted) 0.5" x 9"	Pediatric (Pink) 0.75" x 12"	Adult Reduc. Width (Blue) 0.75" x 18"	Bariatric (Orange) 1.5" x 32"
Smooth	Flat Packed	48NLINB-FP 500/box	48NLINP-FP 500/box	48NLINA-FP 1000/case	48NLPED-FP 1000/case	38TXB-FP 250/pack	*
	Roll Pack	*	*	*	48NLPED-P 25/roll, 40 rolls/case	*	*
	Dispenser Box	*	*	*	48NLPED-DBX 25/roll, 40 rolls/case	38TXB-DBX 25/roll, 40 boxes/case	48NLBAR-DBX 20/roll, 15 rolls/case
	Individually Rolled	*	*	*	48NLPED-IR 480/case	38TXB-IR 1000/case	48NLBAR-O 480/case

*Contact Caresfield to order



8095 215th Street West, Lakeville, Minnesota

952.890.7100 | sales@caresfield.com | caresfield.com

